


**Total Solutions™ for Breast Cancer Diagnosis & Treatment  
 Needles, HDR Afterloader, Applicators...and much more!**



**Best™ Breast Localization Needles**


- Best localization needles are available in 5, 7, 10 and 15 cm lengths
- Available in various ranges of seed activities with a shelf life of 90 days
- Needles are supplied sterile and are individually packaged and shielded

**Best™ HDR Afterloader\***

- Provides 18 channels for dose delivery
- Battery back-up in case of power failure and automated wire recovery
- Source wire tracking using sensors and encoders
- Treatment data maintained in case of failure and treatment resumed after fix
- Quick source replacement process thereby reducing down time
- Available with Cobalt-60 or Iridium-192 sources. Also, Ytterbium-169 source available soon!



**Best™ Double-Balloon Breast Brachytherapy Applicator**



- Improved dose distribution and conformity
- Less dose to critical organs such as skin, lung, heart, chest wall, etc.
- Convenient to use
- US FDA 510(k) registered


\* Product shown is not available for sale currently.

# Total Solutions™ for Breast Cancer Diagnosis & Treatment

## Best™ Sonalis™ Ultrasound Imaging System

The Best™ Sonalis™ Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures.

- Totally sealed, self-healing antimicrobial keyboard with SensoFoil® Technology
- High resolution breast probe provides for needle aspiration, marker insertion and HDR applicator placement
- Longitudinal array on rectal probe provides for 140 mm length of view
- Simultaneous imaging of transverse and sagittal planes
- PC Based System provides a platform for future upgrades and application-specific modules



## Best™ Treatment Planning System (TPS)

Best™ TPS is the planning system of choice for enhanced patient outcomes.

- Can be used for treatment planning with the breast balloon applicator
- Dual Activity Module creates treatment plans using multiple activity sources
- Slice Shifter™ provides ability to correct image slice positioning errors
  - Concurrent 2D & 3D Visualization allows instant visualization of seed placement, dose distribution & anatomical structures
  - Pattern Load Feature – User-definable seed patterns can be used to initiate new plans



## Best™ UCT-1000 Universal Chair


- 4-Function Hand-Held Pendant (Control)
- Electronic variable height control/back support/foot support & horizontal tilt 12 degrees – head or foot
- Padded variable height intravenous arm rest
- Removable Stainless Steel IV Pole
- Stainless steel bilateral cushioned knee supports



## Best™ Brachytherapy Kit

- Catheter tubes are available in 4 colors (purple, green, clear or blue) with either radiopaque or clear nylon buttons
- Half Moon Buttons are available in radiopaque (pictured) or clear nylon

5 Single Leader Catheters



5 Half Moon Buttons



5 Friction Cuffs



5 Red Caps



5 Implant Needles

5 Stylets with Hubs



Flexi Needle Tray

## Best™ Flexi Needles




15G Flexi Needle




16G Flexi Needle



17G Flexi Needle



Stainless Steel or Tungsten Obturator



Flexi Needle with Female Hub



## BATCAM™ Multiprobe Image-Guided Radiation Therapy

BATCAM™ uses the latest compact ultrasound technology to provide increased sensitivity and resolution which yields more clinical information and expanded applications. Plus, it provides the ability to orient the probe in any rotation or translation position.

- Imports treatment plans from all major treatment planning systems
- Enhanced positioning precision using triangulated positioning techniques
- ImageSync™ Interactive Alignment provides improved visualization
- Position verification with image capture ensures accuracy
- User-friendly software components featuring a touch-screen interface








# Best™ medical international

**Total Solutions™ for Breast Cancer Diagnosis & Treatment  
I-125/Pd-103 Seeds, Fiducial Markers...and much more!**


- We manufacture both I-125 & Pd-103 seeds in Springfield, Virginia USA
- Best Medical makes custom seeds and strands to your exact specifications — any source, any configuration, any time!
- Best visualization in ultrasound, fluoro, X-ray, CT and MRI, facilitating real-time dosimetry
- Assistance setting up new programs (licensing, procurement of equipment, proctoring, reimbursement consulting, patient education and marketing)



**Best™ Iodine-125 Seed**



**Best™ Palladium-103 Seed**



**We ship your order within  
24 hours  
Sterile and Non-Sterile**

## Best™ Gold Fiducial Marker Kit



**Stranded**

**Loose**

**Best Medical can  
provide any standard or  
customized gold marker  
configuration made to  
your exact specifications,  
sterile or non-sterile!**

Best Medical International, Inc. © 2014–2018

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA  
tel: 703 451 2378 800 336 4970 [www.besttotalsolutions.com](http://www.besttotalsolutions.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA